# Assays for monitoring susceptibility of influenza viruses to neuraminidase inhibitors

### Margaret Okomo-Adhiambo, Tiffany G. Sheu, Larisa V. Gubareva

Virus Surveillance and Diagnosis Branch, Influenza Division, National Center for Immunization and Respiratory Diseases, Centers for Disease Control and Prevention, Atlanta, GA, USA.

Correspondence: Larisa Gubareva, Centers for Disease Control and Prevention, Mail Stop G-16, 1600 Clifton Road, Atlanta, GA 30333, USA. E-mail: lgubareva@cdc.gov

Close monitoring of drug susceptibility among human influenza viruses was necessitated by widespread resistance to M2 inhibitors in influenza H1N1 (pre-pandemic and 2009 pandemic) and H3N2 viruses, and of oseltamivir resistance in pre-pandemic H1N1 viruses. The FDA-approved neuraminidase (NA) inhibitors (NAIs), oseltamivir and zanamivir, as well as investigational NAIs, peramivir and laninamivir, are currently the principal treatment options for managing influenza infection. However, there are challenges associated with assessing virus susceptibility to this class of drugs. Traditional cell culture-based assays are not reliable for phenotypic testing of NAI susceptibility due to complexity in interpretation. Two types of laboratory assays are currently available for monitoring NAI susceptibility, phenotypic such as the neuraminidase inhibition (NI) assay and genotypic. The NI assay's requirement for propagated virus lengthens testing turnaround; therefore, the need for timely detection of molecular markers associated with NAI resistance (e.g., H275Y in H1N1) has

spurred the development of rapid, high-throughput assays, such as real-time RT-PCR and pyrosequencing. The high sensitivity of genotypic assays allows testing of clinical specimens thus eliminating the need for virus propagation in cell culture. The NI assays are especially valuable when a novel virus emerges or a new NAI becomes available. Modifications continue to be introduced into NI assays, including optimization and data analysis criteria. The optimal assay of choice for monitoring influenza drug susceptibility varies widely depending on the needs of laboratories (e.g., surveillance purposes, clinical settings). Optimally, it is desirable to combine functional and genetic analyses of virus isolates and, when possible, the respective clinical specimens.

**Keywords** Functional analysis, genetic analysis, influenza antiviral susceptibility, neuraminidase inhibitors.

Please cite this paper as: Okomo-Adhiambo et al. (2012) Assays for monitoring susceptibility of influenza viruses to neuraminidase inhibitors. Influenza and Other Respiratory Viruses 7(Suppl. 1), 44–49.

### Introduction

Neuraminidase (NA) inhibitors (NAIs) are currently the only antivirals that are effective for prophylaxis and treatment of influenza A and B infections, due to widespread resistance to M2 inhibitors among currently circulating seasonal influenza A viruses that has greatly diminished the usefulness of this class of drugs. Oseltamivir and zanamivir, administered orally and by inhalation, respectively, are currently the only FDA-approved NAIs for use against type A and type B influenza infections. A newer NAI, peramivir, currently investigational in the United States, has been developed as an intravenous (IV) formulation, and is licensed in Japan And S. Korea, while R-125489, developed as an inhaled prodrug laninamivir (CS-8958), is licensed in Japan.

The 2007–2008 influenza season experienced the dramatic emergence of oseltamivir-resistant seasonal influenza

A(H1N1) viruses carrying the H275Y mutation in the NA, 7-9 which by the 2008–2009 season accounted for up to 100% of seasonal A(H1N1) viruses in many countries. 10,11 These H275Y variants were subsequently replaced by influenza A(H1N1)pdm09 viruses, which emerged in April 2009, leading to a sharp decline in circulating oseltamivir-resistant viruses. In the following 2009-2010 and 2010-2011 seasons, resistance to NAIs among the newly emergent influenza A(H1N1)pdm09 viruses remained low (<1% in the United States). 12,13 However, during the 2010-2011 season, a cluster of oseltamivir-resistant influenza A(H1N1)pdm09 viruses with the H275Y substitution, in circulation between May and September 2011, was reported in Australia. 14,15 The majority of the patients from whom these H275Y variants were recovered had no known oseltamivir exposure. Such potential for emergence and community spread of NAI-resistant variants, coupled with the limited pharmaceutical options currently available for

the control of influenza infections, emphasizes the need to monitor NAI susceptibility among influenza viruses circulating globally.

Assessment of influenza antiviral susceptibility to NAIs is primarily performed using functional phenotypic NA inhibition (NI) assays coupled with genotypic methods such as pyrosequencing 16 and Sanger dideoxy sequence analysis 17 of the NA gene to detect mutations that are associated with NAI resistance. While these methods are effective in detecting influenza viruses resistant to NAIs, there is a need for more high-throughput, affordable and low turnaround approaches for monitoring influenza antiviral susceptibility. This review highlights phenotypic and genotypic approaches that are currently available for assessing influenza virus susceptibility to NAIs.

## Phenotypic methods for assessing influenza virus susceptibility to NAIs

Traditional cell culture–based assays are desirable for initial screening in antiviral susceptibility studies, due to their ability to detect a broad range of resistant phenotypes. However, when applied to influenza, antiviral susceptibility assessed in cell culture lacks correlation with the susceptibility assessed *in vivo* in humans or animal models.<sup>18</sup> In this respect, the NI assay, which functionally assesses the inhibition of the enzyme by the NAI, is beneficial.

Functional methods such as the NI assay allow detection of drug-resistant viruses with established and/or novel changes in the target enzyme. Either the fluorescent<sup>19</sup> or chemiluminescent<sup>20</sup> NI assays are typically the choice for surveillance purposes. Both assays require propagation of virus prior to testing and small synthetic substrates, namely methyl umbelliferone *N*-acetyl neuraminic acid (MUNAN-A)<sup>21</sup> for the fluorescent assay and a 1,2-dioxetane derivative of neuraminic acid<sup>22</sup> for the chemiluminescent assay.

The chemiluminescent and fluorescent NI assays (Table 1) each have advantages and disadvantages associated with their use; for example, the fluorescence-based assay is less costly but requires viruses with higher titers,<sup>23</sup> compared to the chemiluminescence-based assay, which has been shown to provide greater linearity of signal and higher sensitivity in measuring NA activity.<sup>24</sup> The fluorescent assay is preferable for detecting resistance when viral sample permits, as it typically offers better discrimination between NAI susceptible and resistant viruses compared to the chemiluminescent assay.<sup>23</sup> Nevertheless, NAI-resistant mutants can accurately be detected by either NI assay; therefore, the choice of method to use as the primary assay depends on the objectives and requirements of individual surveillance laboratories. Sometimes, an array of assays is applied in characterizing resistance caused by a novel mutation(s).

The NI assay determines the concentration of an NAI needed to reduce enzyme activity by 50% (IC<sub>50</sub>). To determine IC<sub>50</sub> values, raw fluorescent NI assay data (expressed as relative fluorescence units, RFUs) or raw chemiluminescent NI assay data (expressed as relative light units, RLUs) are plotted against drug concentration (nM) using curve fitting software such as JASPR (in-house, CDC)<sup>25</sup> or Robosage (in-house, GlaxoSmithKline). Statistical analyses to determine IC<sub>50</sub> cutoff values for outliers are previously described.  $^{17,25,26}$ 

The reagents used in the chemiluminescent and fluorescent NI assays are commercially available as the NA-Star® Influenza Neuraminidase Inhibitor Reagent Kit and the NA-Fluor<sup>TM</sup> Influenza Neuraminidase Assay Kit, respectively (Applied Biosystems, Foster City, CA, USA). These kits provide validated reagents (except NAIs and reference virus strains) for rapid and sensitive quantitation of influenza NA activity in 96-well microplate formats, enabling improved global assay standardization and more accurate comparison of results between laboratories. The manufacturers' protocols provided in the kits can be optimized to meet individual laboratories' needs. Alternatively, the fluorescent NI assay can be performed using reagents that are prepared in-house using standard chemicals and MUN-ANA (Sigma Aldrich, St. Louis, MO, USA), which is purchased separately. A next-generation chemiluminescencebased assay, the NA-XTD<sup>TM</sup> Influenza Neuraminidase Assay Kit (Applied Biosystems) is also commercially available. This kit provides longer signal readouts compared to the first-generation NA-Star® kit and includes detection reagents that eliminate the need for luminometers equipped with a reagent injector, thereby improving ease-in the-use.

Reference panels of NAI-sensitive and -resistant viruses, to aid in standardizing NI assays and assessing influenza virus susceptibility to NAIs, are available through the ISIRV-Antiviral Group (ISIRV-AVG),<sup>27</sup> the Centers for Disease Control and Prevention (CDC; email: fluantiviral@cdc.gov), or the Influenza Reagent Resource (IRR).

Of note, from a technical standpoint, the NI assay is not a true phenotypic assay and does not account for the interplay of the hemagglutinin (HA) receptor-binding and the NA receptor-destroying activities, which occurs in cell culture. Yet, virus propagation in cell culture is not without limitations, as it may select variants with changes in the virus surface glycoproteins, the HA, and/or the NA, some of which may alter the drug susceptibility profile of the virus in the NI assay. <sup>28,29</sup> Consequently, there is a need for sequence confirmation of both known markers of resistance and changes associated with cell culture selection in addition to functional testing. Nonetheless, virus culture remains an essential component of antiviral resistance monitoring as the NI assay requires virus.

testing
4
=
9
Ę
늄
Ŭ
S
S
=
2
>
Ξ
ä
fluenza antivira
Z
ē
=
ī
-=
5
¥
2
8
$\preceq$
;
ĕ
$\subseteq$
Q
4
0
5
g
ic and genotypic method
$\simeq$
ā
.≍
×
S
n
Ð
유
<del>-</del>
41
-
ā
0

Phenotypic (functional) methods Chemiluminescent NI assay Chemiluminescent NI assay Chemiluminescent NI assay Chemiluminescent NI assay tance Detection Kit NA-Staria Influenza Neuraminidase Assay Kit Assay can be performed using in-house prepared reagents  Genotypic methods Sanger dideoxy sequencing Pyrosequencing Pyrosequencing Sequence analysis (SQA) Shaal-time reverse transcriptase—PCR (RT-PCR) Chemiluminescent NI assay allow accurate detect viruses with established moles wiruses with established moles and viruses with established moles and viruses with assay and viruses of individual surgential surgerial genetic esting allows identify viral genome associated with conferring antiviral resistance Real-time reverse transcriptase—PCR (RT-PCR) Choice of NI assay depends on requirements of individual surgential genotypic testing allows identify conferring antiviral resistance and viral surgential is time-savine.	Advantages  NI assays allow accurate detection of drug-resistant viruses with established molecular markers (e.g., HZ75Y in N1 subtypes) and/or novel changes in the targeted NA enzyme  NI assays provide valuable susceptibility profiles, which cannot be determined solely by genotypic techniques NI assays are available as commercial kits that enable antiviral susceptibility testing to be performed under standardized conditions  Choice of NI assay depends on objectives and requirements of individual surveillance laboratories Genotypic testing can be carried out directly on clinical	Disadvantages/Challenges  Ni testing cannot be carried out directly on clinical material and requires the use of cell grown isolates
nal) methods  NI assay  Iza Neuraminidase Inhibitor Resis-  Iza Neuraminidase Assay Kit  Ara Neuraminidase Assay Kit  formed using in-house  Its  C  C  Lencing  C  C  C  C  C  C  C  C  C  C  C  C  C	curate detection of drug-resistant bilshed molecular markers (e.g., types) and/or novel changes in the rane solely by genotypic techniques able as commercial kits that enable sulity testing to be performed under ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	NI testing cannot be carried out directly on clinical material and requires the use of cell grown isolates
nal) methods  NI assay  Iza Neuraminidase Inhibitor Resis- Tormed using in-house  tits  S (SQA)  polymorphism analysis (SNP)  C anscriptase—PCR (RT-PCR)	urate detection of drug-resistant bilshed molecular markers (e.g., types) and/or novel changes in the me already and/or novel changes in the me already and/or novel changes in the raluable susceptibility profiles, which nined solely by genotypic techniques able as commercial kits that enable illity testing to be performed under ditions depends on objectives and depends on objectives and can be carried out directly on clinical	NI testing cannot be carried out directly on clinical material and requires the use of cell grown isolates
NI assay nza Neuraminidase Inhibitor Resis- nza Neuraminidase Assay Kit nza Neuraminidase Assay Kit formed using in-house nts c c uencing polymorphism analysis (SNP) c c c s (SQA)	vilshed molecular markers (e.g., types) and/or novel changes in the rane and/or novel changes in the raleable susceptibility profiles, which vined solely by genotypic techniques able as commercial kits that enable illity testing to be performed under ditions depends on objectives and andividual surveillance laboratories can be carried out directly on clinical	and requires the use of cell grown isolates
nza Neuraminidase Inhibitor Resis- nza Neuraminidase Assay Kit nza Neuraminidase Assay Kit formed using in-house nts c uencing s (SQA) polymorphism analysis (SNP)	types) and/or novel changes in the rme valuable susceptibility profiles, which valuable solely by genotypic techniques able as commercial kits that enable sility testing to be performed under ditions depends on objectives and andividual surveillance laboratories can be carried out directly on clinical	
nza Neuraminidase Assay Kit  nza Neuraminidase Assay Kit formed using in-house  nts  uencing  s (SQA)  polymorphism analysis (SNP)	rme raluable susceptibility profiles, which ined solely by genotypic techniques sable as commercial kits that enable illity testing to be performed under ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	Elevated IC <sub>50</sub> values must be combined with genotypic
nza Neuraminidase Assay Kit  nza Neuraminidase Assay Kit formed using in-house  nts  uencing  s (SQA)  polymorphism analysis (SNP)	raluable susceptibility profiles, which nined solely by genotypic techniques able as commercial kits that enable illity testing to be performed under ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	information to accurately define resistance
nza Neuraminidase Assay Kit hormed using in-house tts  tts  uencing  polymorphism analysis (SNP)  proscriptase—PCR (RT-PCR)	ined solely by genotypic techniques able as commercial kits that enable illity testing to be performed under ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	There is no established cutoff IC <sub>50</sub> value that is indicative
nza Neuraminidase Assay Kit formed using in-house tts  uencing  s (SQA)  polymorphism analysis (SNP)  snscriptase—PCR (RT-PCR)	able as commercial kits that enable illity testing to be performed under ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	of clinically relevant resistance
formed using in-house  tts  uencing s (SQA) polymorphism analysis (SNP)	uility testing to be performed under ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	Variations in assay conditions may affect IC <sub>50</sub> values
uencing c s (SQA) polymorphism analysis (SNP) c snscriptasa—PCR (RT-PCR)	ditions depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	generated in the NI assay
uencing c (SQA) polymorphism analysis (SNP)	depends on objectives and ndividual surveillance laboratories can be carried out directly on clinical	The fluorescence-based assay requires viruses with higher
uencing C S (SQA) polymorphism analysis (SNP)	ndividual surveillance laboratories can be carried out directly on clinical	titers compared to the chemiluminescence-based assay
uencing  S (SQA)  polymorphism analysis (SNP)  anscriptase—PCR (RT-PCR)	can be carried out directly on clinical	
C rphism analysis (SNP)		Genotypic assays carry a high risk for cross-contamination
C (SQA) tide polymorphism analysis (SNP) transcriptase—PCR (RT-PCR)		Some genotypic methods such as Sanger sequencing are
NP)	Genotypic testing allows identification of mutations on	time-consuming, labor intensive, and expensive
NP)	viral genome associated with amino acid substitutions	Interpretation of mutations is difficult without phenotypic
	al resistance	information
	Genotypic testing is time-saving, high throughput,	Understanding the relationship between phenotype and
coupled with detection methods/chemistries relatively simple, and	relatively simple, and rapid (with exception of Sanger	genotype remains a challenge
such as		Determining the proportion of mixed virus population that
SYBR green agents Methods are not diffic	Methods are not difficult to implement using existing	is defined as resistant is a challenge
		1
Single-nucleotide polymorphism (SNP) analysis Choice of genotypic m	Choice of genotypic method to use depends on	
Hybridization probes	capabilities	
High-resolution melting analysis		
Rolling circle amplification		
Conventional end-point RT-PCR coupled with		
methods such as		
Single-nucleotide polymorphism (SNP)		
genotyping		
Restriction fragment length		
polymorphism (RPLP) analysis		

The NI assay provides valuable quantitative susceptibility data, which cannot be determined solely by sequence-based techniques. The susceptibility of viruses with intermediate IC<sub>50</sub> values is usually difficult to interpret; therefore, such viruses are further investigated to determine the presence of molecular changes in the NA and to determine their frequency of detection as well as potential clinical importance.

The  $IC_{50}$  values generated in NI assays provide valuable information for detecting NAI-resistant viruses, but the lack of an established threshold  $IC_{50}$  value indicative of clinically relevant resistance does not allow  $IC_{50}$ s to be used in drawing direct correlations with drug concentrations required to inhibit virus replication in the infected human host. Nevertheless, the assessment of NAI susceptibility of influenza viruses in the NI assay, reinforced by NA sequence analysis of virus isolates with  $IC_{50}$ s above baseline values and their matching clinical specimens, provides a reliable and reasonably comprehensive approach to the identification of NAI-resistant isolates for surveillance purposes.

# Genotypic methods for assessing influenza virus susceptibility to NAIs

The propensity for rapid and constant evolution of the RNA genome of the influenza virus requires flexible diagnostic tools for monitoring existing and novel drug resistance mutations.

The most commonly applied genotypic methodology for detecting NAI resistance mutations in the NA couples conventional end-point reverse transcription—polymerase chain reaction (RT-PCR) with techniques such as Sanger sequencing and pyrosequencing (Table 1). Full gene sequencing is extremely informative and accurate in detecting changes to the virus genome, but is time-consuming, laborious, and expensive, and may be indiscriminate in determining components of mixed virus populations. Nevertheless, Sanger sequencing remains the assay of choice for identifying both characterized and novel changes, which may underlie phenotypic resistance.

Pyrosequencing technology has for a number of years been applied to the detection of known molecular markers of resistance in the NA gene. 30–33 At the CDC, mutations in the NA that are most commonly screened using pyrosequencing are the oseltamivir resistance conferring H275Y substitution in A(H1N1)pdm09 viruses, the E119V, R292K, and N294S substitutions in influenza A(H3N2) viruses, as well as changes at residues H273, D197, E117, and R374 in influenza B viruses.

Initially used for sequencing target regions (SQA analysis), pyrosequencing assays have been applied in quantification of single-nucleotide polymorphisms (SNP analysis). SNP analysis relies on the premise that most known NAI resistance markers result from single point

mutations in the codons of certain critical residues in the NA. The technology has also been applied in characterizing complex mixtures in the HA.<sup>36</sup>

Pyrosequencing is easily scaled up for high-throughput testing and provides highly informative genetic data for known markers of drug resistance. This approach cannot, however, identify novel changes in the NA that may confer resistance or subtle differences in virus susceptibility to NAIs, and must constantly be updated to accommodate such changes. It is also associated with costly maintenance contracts, specialized equipment, and reagents.

Real-time RT-PCR is a rapid, high-throughput technology commonly used for influenza typing/subtyping in diagnostic and clinical laboratories. Recently, real-time RT-PCR has been applied to drug resistance detection (Table 1), utilizing several detection techniques, such as SYBR green intercalating agents;<sup>37</sup> single-nucleotide polymorphism (SNP) probes;<sup>38–43</sup> high-resolution melting (HRM) analysis;<sup>44,45,46</sup> hybridization probes;<sup>47</sup> minor groove binding (MGB) probes;<sup>48</sup> and rolling circle amplification.<sup>49,50</sup> The ease-of-use, accessibility to equipment, and the availability of reagents and training further contribute to the widespread application of real-time RT-PCR assays.

Additional PCR-based methods, utilizing conventional RT-PCR, have been used for detecting drug resistance mutations in influenza viruses, including SNP analysis by single-nucleotide probe extension (SNaP Shot assay)<sup>35</sup> and restriction fragment length polymorphism (RFLP) analysis.<sup>51,52</sup>

Regardless of the RT-PCR method or reaction chemistry, benefits of PCR-based characterization include lower cost, potential implementation in a high-throughput system, and generally straightforward interpretation.

The main risk of genotypic tests is cross-contamination, requiring specific procedures for its prevention. The risk of false-negative results due to insufficient extraction or procedural errors still remains. Less frequent causes of genotypic test error are mixtures (wild-type strain plus emerging mutant or multiple mutations) and silent mutations, with change of nucleotide but not amino acid.

#### Conclusion

The choice of assay for assessing influenza virus susceptibility to NAIs depends on factors pertaining to appropriateness to the setting, cost, sustainability, speed in obtaining valid results, reliability in terms of predictive values, and accessibility. Although slow conventional tests could still be used for epidemiological monitoring of drug resistance, rapid genotypic testing facilitates more appropriate patient management and can significantly advance large-scale epidemiological studies of drug-resistant variants. In practice, it is likely that more than one method will be needed

to ensure rapid as well as accurate detection of resistance to NAIs.

### **Conflicts of interest**

The findings and conclusions of this report are those of the authors and do not necessarily represent the views of the Centers for Disease Control and Prevention (CDC). The authors have no potential conflicts to declare.

#### References

- **1** Bright RA, Shay DK, Shu B, Cox NJ, Klimov AI. Adamantane resistance among influenza A viruses isolated early during the 2005–2006 influenza season in the United States. JAMA 2006; 295:891–894.
- 2 Deyde VM, Xu X, Bright RA et al. Surveillance of resistance to adamantanes among influenza A(H3N2) and A(H1N1) viruses isolated worldwide. J Infect Dis 2007; 196:249–257.
- **3** Gubareva LV, Trujillo AA, Okomo-Adhiambo M *et al.* Comprehensive assessment of 2009 pandemic influenza A (H1N1) virus drug susceptibility in vitro. Antivir Ther 2010; 15:1151–1159.
- **4** Shetty AK, Peek LA. Peramivir for the treatment of influenza. Expert Rev Anti Infect Ther 2012; 10:123–143.
- 5 Sugaya N, Kohno S, Ishibashi T, Wajima T, Takahashi T. Efficacy, safety, and pharmacokinetics of intravenous peramivir in children with 2009 pandemic H1N1 influenza A virus infection. Antimicrob Agents Chemother 2012; 56:369–377.
- 6 Sunagawa S, Higa F, Cash HL, Tateyama M, Uno T, Fujita J. Single-dose inhaled laninamivir: registered in Japan and its potential role in control of influenza epidemics. Influenza Other Respi Viruses 2012; DOI: 10.1111/j.1750-2659.2012.00351.x..
- 7 CDC. Influenza activity United States and worldwide, May 18-September 19, 2008. MMWR Morb Mortal Wkly Rep 2008; 57:1046–1049.
- 8 Lackenby A, Hungnes O, Dudman SG et al. Emergence of resistance to oseltamivir among influenza A(H1N1) viruses in Europe. Euro Surveill 2008; 13(5):8026.
- **9** Hauge SH, Dudman S, Borgen K, Lackenby A, Hungnes O. Oseltamivir-resistant influenza viruses A (H1N1), Norway, 2007–08. Emerg Infect Dis 2009; 15:155–162.
- 10 CDC. Update: influenza activity United States, September 28, 2008–April 4, 2009, and composition of the 2009–10 influenza vaccine. MMWR Morb Mortal Wkly Rep 2009; 58:369–374.
- 11 Baranovich T, Saito R, Suzuki Y et al. Emergence of H274Y oseltamivir-resistant A(H1N1) influenza viruses in Japan during the 2008– 2009 season. J Clin Virol 2010; 47:23–28.
- **12** Graitcer SB, Gubareva L, Kamimoto L *et al.* Characteristics of patients with oseltamivir-resistant pandemic (H1N1) 2009, United States. Emerg Infect Dis 2011; 17:255–257.
- 13 Storms AD, Gubareva LV, Su S et al. Oseltamivir-resistant pandemic (H1N1) 2009 virus infections, United States, 2010–11. Emerg Infect Dis 2012; 18:308–311.
- 14 Hurt AC, Hardie K, Wilson NJ et al. Community transmission of oseltamivir-resistant A(H1N1)pdm09 influenza. N Engl J Med 2011; 365:2541–2542.
- 15 Hurt AC, Hardie K, Wilson NJ et al. Characteristics of a widespread community cluster of H275Y oseltamivir-resistant A(H1N1)pdm09 influenza in Australia. J Infect Dis 2012; 206:148–157. Epub 2012 May 4.

- 16 Deyde VM, Gubareva LV. Influenza genome analysis using pyrosequencing method: current applications for a moving target. Expert Rev Mol Diagn 2009; 9:493–509.
- **17** Sheu TG, Deyde VM, Okomo-Adhiambo M *et al.* Surveillance for neuraminidase inhibitor resistance among human influenza A and B viruses circulating worldwide from 2004 to 2008. Antimicrob Agents Chemother 2008; 52:3284–3292.
- 18 Tisdale M. Monitoring of viral susceptibility: new challenges with the development of influenza NA inhibitors. Rev Med Virol 2000; 10:45–55.
- **19** Hurt AC, Okomo-Adhiambo M, Gubareva LV. The fluorescence neuraminidase inhibition assay: a functional method for detection of influenza virus resistance to the neuraminidase inhibitors. Methods Mol Biol 2012; 865:115–125.
- 20 Okomo-Adhiambo M, Hurt AC, Gubareva LV. The chemiluminescent neuraminidase inhibition assay: a functional method for detection of influenza virus resistance to the neuraminidase inhibitors. Methods Mol Biol 2012; 865:95–113.
- 21 Potier M, Mameli L, Belisle M, Dallaire L, Melancon SB. Fluorometric assay of neuraminidase with a sodium (4-methylumbelliferyl-alpha-D-N-acetylneuraminate) substrate. Anal Biochem 1979; 94:287–296.
- 22 Buxton RC, Edwards B, Juo RR, Voyta JC, Tisdale M, Bethell RC. Development of a sensitive chemiluminescent neuraminidase assay for the determination of influenza virus susceptibility to zanamivir. Anal Biochem 2000; 280:291–300.
- 23 Nguyen HT, Sheu TG, Mishin VP, Klimov AI, Gubareva LV. Assessment of pandemic and seasonal influenza A (H1N1) virus susceptibility to neuraminidase inhibitors in three enzyme activity inhibition assays. Antimicrob Agents Chemother 2010; 54:3671–3677.
- 24 Wetherall NT, Trivedi T, Zeller J et al. Evaluation of neuraminidase enzyme assays using different substrates to measure susceptibility of influenza virus clinical isolates to neuraminidase inhibitors: report of the neuraminidase inhibitor susceptibility network. J Clin Microbiol 2003; 41:742–750.
- 25 Okomo-Adhiambo M, Sleeman K, Ballenger K et al. Neuraminidase inhibitor susceptibility testing in human influenza viruses: a laboratory surveillance perspective. Viruses 2010; 2:2269–2289.
- 26 Monto AS, McKimm-Breschkin JL, Macken C et al. Detection of influenza viruses resistant to neuraminidase inhibitors in global surveillance during the first 3 years of their use. Antimicrob Agents Chemother 2006; 50:2395–2402.
- 27 ISIRV-AVG. Available at http://www.isirv.org/site/index.php/reference-panel (Accessed 31 October 2012).
- 28 Hurt AC, Holien JK, Parker M, Kelso A, Barr IG. Zanamivir-resistant influenza viruses with a novel neuraminidase mutation. J Virol 2009; 83:10366–10373.
- **29** Okomo-Adhiambo M, Nguyen HT, Sleeman K *et al.* Host cell selection of influenza neuraminidase variants: implications for drug resistance monitoring in A(H1N1) viruses. Antiviral Res 2010; 85:381–388
- 30 Deyde VM, Okomo-Adhiambo M, Sheu TG et al. Pyrosequencing as a tool to detect molecular markers of resistance to neuraminidase inhibitors in seasonal influenza A viruses. Antiviral Res 2009; 81:16–24.
- **31** Deyde VM, Sheu TG, Trujillo AA *et al.* Detection of molecular markers of drug resistance in 2009 pandemic influenza A (H1N1) viruses by pyrosequencing. Antimicrob Agents Chemother 2010; 54:1102–1110.
- **32** Sheu TG, Deyde VM, Garten RJ, Klimov Al, Gubareva LV. Detection of antiviral resistance and genetic lineage markers in influenza B virus neuraminidase using pyrosequencing. Antiviral Res 2010; 85:354–360.

- 33 Bao JR, Huard TK, Piscitelli AE et al. Reverse-transcription polymerase chain reaction/pyrosequencing to characterize neuraminidase H275 residue of influenza A 2009 H1N1 virus for rapid and specific detection of the viral oseltamivir resistance marker in a clinical laboratory. Diagn Microbiol Infect Dis 2011; 71:396– 402
- **34** Lackenby A, Democratis J, Siqueira MM, Zambon MC. Rapid quantitation of neuraminidase inhibitor drug resistance in influenza virus quasispecies. Antivir Ther 2008; 13:809–820.
- **35** Duan S, Boltz DA, Li J *et al.* Novel genotyping and quantitative analysis of neuraminidase inhibitor resistance-associated mutations in influenza a viruses by single-nucleotide polymorphism analysis. Antimicrob Agents Chemother 2011; 55:4718–4727.
- 36 Levine M, Sheu TG, Gubareva LV, Mishin VP. Detection of hemagglutinin variants of the pandemic influenza A (H1N1) 2009 virus by pyrosequencing. J Clin Microbiol 2011; 49:1307–1312.
- 37 Medina RA, Rojas M, Tuin A et al. Development and characterization of a highly specific and sensitive SYBR green reverse transcriptase PCR assay for detection of the 2009 pandemic H1N1 influenza virus on the basis of sequence signatures. J Clin Microbiol 2011; 49:335–344
- **38** Wong S, Pabbaraju K, Wong A, Fonseca K, Drews SJ. Development of a real-time RT-PCR assay for detection of resistance to oseltamivir in influenza A pandemic (H1N1) 2009 virus using single nucleotide polymorphism probes. J Virol Methods 2011; 173:259–265.
- **39** Operario DJ, Moser MJ, St GK. Highly sensitive and quantitative detection of the H274Y oseltamivir resistance mutation in seasonal A/H1N1 influenza virus. J Clin Microbiol 2010; 48:3517–3524.
- 40 Renaud C, Kuypers J, Corey L. Diagnostic accuracy of an allele-specific reverse transcriptase-PCR assay targeting the H275Y oseltamivir resistant mutation in 2009 pandemic influenza A/H1N1 virus. J Clin Virol 2010: 49:21–25.
- **41** van der Vries E, Jonges M, Herfst S *et al.* Evaluation of a rapid molecular algorithm for detection of pandemic influenza A (H1N1) 2009 virus and screening for a key oseltamivir resistance (H275Y) substitution in neuraminidase. J Clin Virol 2010; 47:34–37.
- 42 Bolotin S, Robertson AV, Eshaghi A et al. Development of a novel real-time reverse-transcriptase PCR method for the detection of H275Y positive influenza A H1N1 isolates. J Virol Methods 2009; 158:190–194.

- **43** Carr MJ, Sayre N, Duffy M, Connell J, Hall WW. Rapid molecular detection of the H275Y oseltamivir resistance gene mutation in circulating influenza A (H1N1) viruses. J Virol Methods 2008; 153:257–262.
- 44 Tong SY, Dakh F, Hurt AC et al. Rapid detection of the H275Y osel-tamivir resistance mutation in influenza A/H1N1 2009 by single base pair RT-PCR and high-resolution melting. PLoS ONE 2011; 6:e21446
- **45** Lee HK, Lee CK, Loh TP, Tang JW, Tambyah PA, Koay ES. High-resolution melting approach to efficient identification and quantification of H275Y mutant influenza H1N1/2009 virus in mixed-virus-population samples. J Clin Microbiol 2011; 49:3555–3559.
- 46 Redlberger-Fritz M, Aberle SW, Strassl R, Popow-Kraupp T. Rapid identification of neuraminidase inhibitor resistance mutations in seasonal influenza virus A(H1N1), A(H1N1)2009, and A(H3N2) subtypes by melting point analysis. Eur J Clin Microbiol Infect Dis 2012; 31:1593–1601.
- **47** Whiley DM, Jacob K, Nakos J *et al.* Improved detection of genetic markers of antimicrobial resistance by hybridization probe-based melting curve analysis using primers to mask proximal mutations: examples include the influenza H275Y substitution. J Antimicrob Chemother 2012; 67:1375–1379.
- **48** Hindiyeh M, Ram D, Mandelboim M *et al.* Rapid detection of influenza A pandemic (H1N1) 2009 virus neuraminidase resistance mutation H275Y by real-time reverse transcriptase PCR. J Clin Microbiol 2010; 48:1884–1887.
- **49** Wang B, Dwyer DE, Blyth CC *et al.* Detection of the rapid emergence of the H275Y mutation associated with oseltamivir resistance in severe pandemic influenza virus A/H1N1 09 infections. Antiviral Res 2010: 87:16–21.
- **50** Steain MC, Dwyer DE, Hurt AC *et al.* Detection of influenza A H1N1 and H3N2 mutations conferring resistance to oseltamivir using rolling circle amplification. Antiviral Res 2009; 84:242–248.
- **51** Guo L, Garten RJ, Foust AS *et al.* Rapid identification of oseltamivir-resistant influenza A(H1N1) viruses with H274Y mutation by RT-PCR/restriction fragment length polymorphism assay. Antiviral Res 2009; 82:29–33.
- **52** Nukiwa N, Suzuki A, Furuse Y *et al.* Simplified screening method for detecting oseltamivir resistant pandemic influenza A (H1N1) 2009 virus by a RT-PCR/restriction fragment length polymorphism assay. J Virol Methods 2010; 170:165–168.